

Qi400

30Litre Ultrasonic Bath

Ultrawave's **Qi-Series** are robust and reliable, delivering outstanding cleaning performance in a wide range of industries. From precision, pharmaceutical and medical device manufacturing to heavy duty cleaning in automotive and engineering processing. **Qi-Series** Ultrasonic Cleaners are Class I Medical Devices, compliant to MDR (EU) 2017/745.

The 2mm thick stainless steel tank provides longevity and durability whilst excellent ultrasonic activity is generated via Ultrawave's *Frequency Leap* technology.

Frequency Leap technology reduces dead spots and promotes homogenous cavitation activity within the tank. This provides a more consistent scrubbing action for highly effective cleaning results.



The **Qi-400** is a 30Litre ultrasonic cleaner with a menu driven, LCD control panel. This allows simple adjustment of cycle time (0 to 99mins or infinity mode), temperature (ambient to 80°C), ultrasonic power and frequency mode for precise and accurate processing.

Incorporating features such as degas mode and a fluid level sensor, the

Qi-Series is reliable and durable, giving you consistent cleaning results every time.



Specifications are subject to change for continual improvements

Qi400

30Litre Ultrasonic Bath

| | |
|--|---|
| Construction Case | 304 Stainless Steel |
| Construction Tank | 316 Stainless Steel |
| Control | Digital LCD |
| Overall Dimensions | 530 x 380 x 395mm (WxDxH) |
| Tank Dimensions | 485 x 275 x 220mm (WxDxH) |
| Internal Basket Dimensions | 455 x 250 x 140mm (WxDxH) |
| Maximum Capacity | 30Litres |
| Working Capacity | 26Litres |
| Ultrasonic Power | 400W |
| Ultrasonic Power (per litre) | 15W/Litre |
| Operating Frequency | 32 - 38kHz <i>Frequency Leap</i> |
| Maximum Heating Capacity | Ambient - 80°C |
| Heater Power | 1350W |
| Timer | 0 - 99 minutes or Continuous |
| Mains Connection | 230VAC |
| Drain Outlet | Lockable 1/2" BSP valve (located at rear of machine) |
| Lid | 304 Stainless Steel |
| Basket | 304 Stainless Steel |
| Gross Weight (inc. lid, basket & packaging) | 24.5kg |
| Net Weight | 22kg |
| Warranty | 1 Year as standard - extend to 3 Years if registered online with us (within 30days of purchase) |